

# Rotating Disc Electrode

## Work-Up and Questions

- 1) Determine the limiting current as a function of rotation rate and construct a graph of  $\omega^{-1/2}$  versus  $i^{-1}$ .

Figure 1 shows a plot of current versus time, and the Koutecky-Levitch plot obtained from the data, remembering that  $\omega=2\pi f$ .




Figure 1. (left) Data at varied rotation rates using the two solutions. (right) Koutecky-Levitch plot obtained from the data.

- 2) Use linear regression to determine the values for  $i_k$  and  $D_o$ , assuming  $\nu = 1.0074 \times 10^{-6} \text{ m}^2 \text{ s}^{-1}$  for these solutions.

The slope of the line in the Koutecky-Levitch plot was  $466 \text{ s}^{1/2} \text{ A}^{-1}$ . The intercept was  $-7.4 \times 10^{-1} \text{ A}^{-1}$ . The intercept yields a value of  $-1.3 \text{ A}$  for  $i_k$ . The electrode was a Pt disk with diameter about 5 mm.

To calculate  $D_o$ ,

$$\text{slope} = \frac{1}{0.62nFAD_o^{2/3}\nu^{-1/6}C_o^*}$$

$$D_o^{2/3} = \frac{\nu^{1/6}}{\text{slope} \times 0.62nFAC_o^*}$$

$$D_o^{2/3} = \frac{(1.0074 \times 10^{-6})^{1/6} \text{m}^{1/3} \times \frac{100^{1/3} \text{cm}^{1/3}}{\text{m}^{1/3}}}{\frac{466 \text{s}^{1/2}}{\text{A}} \times \text{s}^{1/6} \times \frac{1 \text{As}}{\text{C}} \times 0.62 \times 1 \times \frac{96485 \text{C}}{\text{mol}} \times 0.196 \text{cm}^2 \times \frac{0.15 \text{mol}}{1000 \text{cm}^3}}$$

$$D_o^{2/3} = \frac{0.100 \times 4.64 \text{cm}^{4/3}}{819.6 \text{s}^{3/2}}$$

$$D_O = 1.4 \times 10^{-5} \text{ cm}^2 \text{ s}^{-1}$$

The value according to the Handbook of Chemistry and Physics is  $1.18 \times 10^{-5} \text{ cm}^2 \text{ s}^{-1}$ , quite close to that obtained from these measurements.