

# Cyclic Voltammetry

## Chemicals and Instruments

N<sub>2</sub>-purged solution of 0.5 M H<sub>2</sub>SO<sub>4</sub>

N<sub>2</sub>-purged 2 mM solution of *N,N,N',N'*-tetramethyl-*p*-phenylene diamine in 0.5 M H<sub>2</sub>SO<sub>4</sub>

Potentiostat

2 Pt electrodes

Ag/AgCl/KCl(sat'd) reference electrode

Electrochemical cell or 3-necked flask


## Electrode Reactions with Coupled Homogeneous Chemical Reactions (B&F Chapters 6 and 12, also Nicholson & Shain)

### Task

Use cyclic voltammetry to derive a reaction pathway for the electrooxidation of *N,N,N',N'*-tetramethyl-*p*-phenylene diamine.

### Fundamentals

In homework set 3 (B&F question 6.9), we analyzed cyclic voltammetry data for azotoluene. In this laboratory exercise we will study cyclic voltammetry for *N,N,N',N'*-tetramethyl-*p*-phenylene diamine and apply the analysis techniques we have learned using our experimental data. Figure 1 shows the structure of the molecule and examples of typical cyclic voltammograms with the upper potential limit varied.



**Figure 1.** (left) Structure of *N,N,N',N'*-tetramethyl-*p*-phenylene diamine. (center and right) Cyclic voltammograms of a 2 mM solution of *N,N,N',N'*-tetramethyl-*p*-phenylene diamine in 0.5 M H<sub>2</sub>SO<sub>4</sub> (solid curves) with CVs for the electrolyte solution shown with dashed lines. The scan rate was 100 mV/s for both.

## Procedure

Use a Pt working electrode, Pt counter electrode, a Ag/AgCl/KCl (sat'd) reference electrode and a potentiostat to record cyclic voltammograms for the solution of *N,N,N',N'*-tetramethyl-*p*-phenylene diamine at varied scan rates and for the two scan windows shown in the Figure 1.

Try to get your range of scan rates to cover at least a factor of 10, as Nicholson & Shain and Prof. Bocarsly tell us to do. Since the maximum scan window looks to be 0–1.6 V versus Ag/AgCl, a complete scan at 20 mV/s will require less than 3 minutes, so that seems like a reasonable lower limit for the scan rate (provided data also looks good at 200 mV/s).

Collect data using at least 5 scan rates, the more data you collect the better you will be able to fit the current functions.

### Work-Up and Questions

- 1) Plot your cyclic voltammetry data for each scan rate, including overlays for the baseline scans. To do this neatly you may need a separate plot for each scan-window/scan-rate combination.
- 2) Use your data to make a table including  $i_{pc}$ ,  $i_{pa}$ ,  $E_{pc}$ ,  $E_{pa}$ , and  $E_{p/2}$  (as applicable), similar to the table that accompanied B&F question 6.9.
- 3) From the peak ratios and dependence of peak position on scan rate, comment on whether each wave appears reversible or irreversible.
- 4) Use appropriate plots and fits to analyze the data. Assume  $n=1$ . For each reversible wave, determine a diffusion coefficient. For each irreversible wave, determine  $\alpha$ , and a diffusion coefficient.
- 5) For each wave, which mechanism from Nicholson & Shain best fits the data?
- 6) Speculate as to a possible chemical mechanism for the observed behavior.

Adapted from Holze, R. *Experimental Electrochemistry A Laboratory Textbook*. Wiley-VCH, Darmstadt. 2009.